



Product Datasheet

CD41 antibody (orb43845)

Description	Mouse monoclonal antibody to CD41
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr
Target	CD41b
Preservatives	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG3
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P08514
Entrez	3674
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of human CD42b positive debris...



Flow cytometry surface staining pattern ...