

Product Datasheet

Human CCL14 ELISA Kit (orb438450)

Description	Human Chemokine C-C-Motif Ligand 14 (CCL14) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q16627
Sensitivity	6.8pg/mL
Entrez	6358
Expiration Date	Please enquire.