

Product Datasheet

Human CCL1 ELISA Kit (orb438451)

Description	Human Chemokine C-C-Motif Ligand 1 (CCL1) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P22362
Sensitivity	27pg/mL
Entrez	6346
Expiration Date	Please enquire.