

Product Datasheet

Mouse CCL1 ELISA Kit (orb438452)

Description	Mouse Chemokine C-C-Motif Ligand 1 (CCL1) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P10146
Sensitivity	32pg/mL
Entrez	20290
Expiration Date	Please enquire.