

Product Datasheet

Human CCNB2 ELISA Kit (orb438456)

Description	Human Cyclin B2 (CCNB2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O95067
Sensitivity	0.066ng/mL
Entrez	9133
Expiration Date	Please enquire.