

Product Datasheet

Mouse CCND1 ELISA Kit (orb438459)

Description	Mouse Cyclin D1 (CCND1) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P25322
Sensitivity	5.7pg/mL
Entrez	12443
Expiration Date	Please enquire.