

Product Datasheet

Human CCND2 ELISA Kit (orb438461)

Description	Human Cyclin D2 (CCND2) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P30279
Sensitivity	5.7pg/mL
Entrez	894
Expiration Date	Please enquire.