

Product Datasheet

Mouse CCR10 ELISA Kit (orb438462)

Description	Mouse Chemokine C-C-Motif Receptor 10 (CCR10) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9JL21
Sensitivity	0.051ng/mL
Entrez	12777
Expiration Date	Please enquire.