



Product Datasheet

CD235a antibody (orb43847)

Description	Mouse Monoclonal to CD235a.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	AGG, FC, IHC-Fr, IHC-P
Immunogen	Synthetic peptide (Human, N-terminal)
Target	CD235ab
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG2b
Clonality	Monoclonal
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Flow cytometry analysis (surface staining)



Immunohistochemistry staining of human aorta