

Product Datasheet

Mouse CCT2 ELISA Kit (orb438474)

Description	Mouse Chaperonin Containing TCP1, Subunit 2 (CCT2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P80314
Sensitivity	0.058ng/mL
Entrez	12461
Expiration Date	Please enquire.