

Product Datasheet

Human CD14 ELISA Kit (orb438476)

Description	Human Cluster Of Differentiation 14 (CD14) Sandwich ELISA Kit
Reactivity	Human
Range	9.375-600ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, cell lysates and other biological fluids.
Sample Types	serum, plasma, cell lysates and other biological fluids.
Uniprot ID	P08571
Sensitivity	3.3ng/mL
Entrez	929
Expiration Date	Please enquire.