

Product Datasheet

Mouse CD14 ELISA Kit (orb438477)

Description	Mouse Cluster Of Differentiation 14 (CD14) Sandwich ELISA Kit
Reactivity	Mouse
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P10810
Sensitivity	0.72ng/mL
Entrez	12475
Expiration Date	Please enquire.