

Product Datasheet

Human CD160 ELISA Kit (orb438480)

Description	Human Cluster Of Differentiation 160 (CD160) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, cell lysates and other biological fluids.
Sample Types	serum, plasma, cell lysates and other biological fluids.
Uniprot ID	O95971
Sensitivity	30pg/mL
Entrez	11126
Expiration Date	Please enquire.