

Product Datasheet

Rat CD163 ELISA Kit (orb438482)

Description	Rat Cluster Of Differentiation (CD163) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	D3Z9U2
Sensitivity	0.113ng/mL
Entrez	312701
Expiration Date	Please enquire.