

Product Datasheet

Mouse CD19 ELISA Kit (orb438484)

Description	Mouse Cluster Of Differentiation 19 (CD19) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P25918
Sensitivity	0.047ng/mL
Entrez	12478
Expiration Date	Please enquire.