


Product Datasheet

CD235a antibody (PE) (orb43849)

Description	Mouse Monoclonal to CD235a.
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	Synthetic peptide (Human, N-terminal)
Target	CD235ab
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 20 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.
Isotype	Mouse IgG2b
Clonality	Monoclonal
Dilution Range	Flow cytometry: The reagent is designed for analysis of human blood cells using 20 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...