

Product Datasheet

Human CD27BP ELISA Kit (orb438492)

Description	Human CD27 Binding Protein (CD27BP) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O15304
Sensitivity	0.068ng/mL
Entrez	10572
Expiration Date	Please enquire.