

Product Datasheet

Chicken CD30L ELISA Kit (orb438493)

Description	Chicken Cluster Of Differentiation 30 Ligand (CD30L) Sandwich ELISA Kit
Reactivity	Gallus
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q800J1
Sensitivity	0.034ng/mL
Entrez	378930
Expiration Date	Please enquire.