

Product Datasheet

Human CD30L ELISA Kit (orb438494)

Description	Human Cluster Of Differentiation 30 Ligand (CD30L) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P32971
Sensitivity	33pg/mL
Entrez	944
Expiration Date	Please enquire.