

Product Datasheet

Mouse CD34 ELISA Kit (orb438497)

Description	Mouse Cluster Of Differentiation 34 (CD34) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q64314
Sensitivity	0.061ng/mL
Entrez	12490
Expiration Date	Please enquire.