

Product Datasheet

Bovine CD40L ELISA Kit (orb438499)

Description	Bovine Cluster Of Differentiation 40 Ligand (CD40L) Sandwich ELISA Kit
Reactivity	Bovine
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P51749
Sensitivity	0.18ng/mL
Entrez	282387
Expiration Date	Please enquire.