

Product Datasheet

Canine CD40L ELISA Kit (orb438500)

Description	Canine Cluster Of Differentiation 40 Ligand (CD40L) Sandwich ELISA Kit
Reactivity	Canine
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O97626
Sensitivity	0.032ng/mL
Entrez	403468
Expiration Date	Please enquire.