

Product Datasheet

Porcine CD40L ELISA Kit (orb438503)

Description	Porcine Cluster Of Differentiation 40 Ligand (CD40L) Sandwich ELISA Kit
Reactivity	Porcine
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q95MQ5
Sensitivity	0.054ng/mL
Entrez	397231
Expiration Date	Please enquire.