

Product Datasheet

Rat CD40L ELISA Kit (orb438504)

Description	Rat Cluster Of Differentiation 40 Ligand (CD40L) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9Z2V2
Sensitivity	0.061ng/mL
Entrez	84349
Expiration Date	Please enquire.