

Product Datasheet

Human CD59 ELISA Kit (orb438507)

Description	Human Protectin (CD59) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, cell lysates and other biological fluids.
Sample Types	serum, plasma, cell lysates and other biological fluids.
Uniprot ID	P13987
Sensitivity	13.7pg/mL
Entrez	966
Expiration Date	Please enquire.