

Product Datasheet

Human CD5L ELISA Kit (orb438509)

Description	Human CD5 Antigen Like Protein (CD5L) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O43866
Sensitivity	28.6pg/mL
Entrez	922
Expiration Date	Please enquire.