

Product Datasheet

Mouse CD5L ELISA Kit (orb438510)

Description	Mouse CD5 Antigen Like Protein (CD5L) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9QWK4
Sensitivity	0.136ng/mL
Entrez	11801
Expiration Date	Please enquire.