

Product Datasheet

Human CD83 ELISA Kit (orb438512)

Description	Human Cluster Of Differentiation 83 (CD83) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, cell lysates and other biological fluids.
Sample Types	serum, plasma, cell lysates and other biological fluids.
Uniprot ID	Q01151
Sensitivity	13.2pg/mL
Entrez	9308
Expiration Date	Please enquire.