

Product Datasheet

Mouse CD83 ELISA Kit (orb438513)

Description	Mouse Cluster Of Differentiation 83 (CD83) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O88324
Sensitivity	35pg/mL
Entrez	12522
Expiration Date	Please enquire.