

Product Datasheet

Human CDCP1 ELISA Kit (orb438520)

Description	Human CUB Domain Containing Protein 1 (CDCP1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9H5V8
Sensitivity	0.127ng/mL
Entrez	64866
Expiration Date	Please enquire.