

Product Datasheet

Human CDH5 ELISA Kit (orb438526)

Description	Human Cadherin 5 (CDH5) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P33151
Sensitivity	0.52ng/mL
Entrez	1003
Expiration Date	Please enquire.