

Product Datasheet

Bovine CDHE ELISA Kit (orb438529)

Description	Bovine Cadherin, Epithelial (CDHE) Sandwich ELISA Kit
Reactivity	Bovine
Range	0.937-60ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q6R8F2
Sensitivity	0.23ng/mL
Entrez	282637
Expiration Date	Please enquire.