

Product Datasheet

Porcine CDHE ELISA Kit (orb438532)

Description	Porcine Cadherin, Epithelial (CDHE) Sandwich ELISA Kit
Reactivity	Porcine
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	C6EVT4
Sensitivity	0.195ng/mL
Entrez	100048953
Expiration Date	Please enquire.