

Product Datasheet

Human CDK16 ELISA Kit (orb438538)

Description	Human Cyclin Dependent Kinase 16 (CDK16) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q00536
Sensitivity	12.6pg/mL
Entrez	5127
Expiration Date	Please enquire.