

Product Datasheet

Rat CDKN1A ELISA Kit (orb438549)

Description	Rat Cyclin Dependent Kinase Inhibitor 1A (CDKN1A) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q64315
Sensitivity	0.068ng/mL
Entrez	114851
Expiration Date	Please enquire.