

Product Datasheet

Human CEACAM6 ELISA Kit (orb438558)

Description	Human Carcinoembryonic Antigen Related Cell Adhesion Molecule 6 (CEACAM6) Sandwich ELISA Kit
Reactivity	Human
Range	23.437-1500pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P40199
Sensitivity	9.68pg/mL
Entrez	4680
Expiration Date	Please enquire.