

Product Datasheet

Human CEBPa ELISA Kit (orb438564)

Description	Human CCAAT/Enhancer Binding Protein Alpha (CEBPa) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P49715
Sensitivity	0.055ng/mL
Entrez	1050
Expiration Date	Please enquire.