

Product Datasheet

Mouse CEBPb ELISA Kit (orb438566)

Description	Mouse CCAAT/Enhancer Binding Protein Beta (CEBPb) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P28033
Sensitivity	33pg/mL
Entrez	12608
Expiration Date	Please enquire.