

Product Datasheet

Human CECR1 ELISA Kit (orb438567)

Description	Human Cat Eye Syndrome Chromosome Region, Candidate 1 (CECR1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9NZK5
Sensitivity	0.115ng/mL
Entrez	51816
Expiration Date	Please enquire.