

Product Datasheet

Human CELSR2 ELISA Kit (orb438568)

Description	Human Cadherin EGF LAG Seven Pass G-Type Receptor 2 (CELSR2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9HCU4
Sensitivity	0.064ng/mL
Entrez	1952
Expiration Date	Please enquire.