

Product Datasheet

Rat CELSR2 ELISA Kit (orb438569)

Description	Rat Cadherin EGF LAG Seven Pass G-Type Receptor 2 (CELSR2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9QYP2
Sensitivity	0.24ng/mL
Entrez	83465
Expiration Date	Please enquire.