

Product Datasheet

Human CEMP1 ELISA Kit (orb438570)

Description	Human Cementum Protein 1 (CEMP1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q6PRD7
Sensitivity	0.124ng/mL
Entrez	752014
Expiration Date	Please enquire.