

Product Datasheet

Human CER ELISA Kit (orb438573)

Description	Human Cerberus (CER) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	095813
Sensitivity	12.3pg/mL
Entrez	9350
Expiration Date	Please enquire.