

Product Datasheet

Rat CES1 ELISA Kit (orb438575)

Description	Rat Carboxylesterase 1 (CES1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q63108
Sensitivity	0.107ng/mL
Entrez	29225
Expiration Date	Please enquire.