

Product Datasheet

Rabbit CETP ELISA Kit (orb438577)

Description	Rabbit Cholesteryl Ester Transfer Protein (CETP) Sandwich ELISA Kit
Reactivity	Rabbit
Range	93.75-6000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P22687
Sensitivity	42.61pg/mL
Entrez	100327267
Expiration Date	Please enquire.