

Product Datasheet

Rat CF6 ELISA Kit (orb438579)

Description	Rat Coupling Factor 6 (CF6) Sandwich ELISA Kit
Reactivity	Rat
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P21571
Sensitivity	0.29ng/mL
Entrez	94271
Expiration Date	Please enquire.