

Product Datasheet

Rat CFB ELISA Kit (orb438585)

Description	Rat Complement Factor B (CFB) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q6MG74
Sensitivity	8.8ng/mL
Entrez	294257
Expiration Date	Please enquire.