

Product Datasheet

Mouse CFH ELISA Kit (orb438590)

Description	Mouse Complement Factor H (CFH) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P06909
Sensitivity	35pg/mL
Entrez	12628
Expiration Date	Please enquire.