

Product Datasheet

Human CFLAR ELISA Kit (orb438599)

Description	Human CASP8 And FADD Like Apoptosis Regulator (CFLAR) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O15519
Sensitivity	0.061ng/mL
Entrez	8837
Expiration Date	Please enquire.