

Product Datasheet

Mouse CHAF1B ELISA Kit (orb438609)

Description	Mouse Chromatin Assembly Factor 1, Subunit B (CHAF1B) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9D0N7
Sensitivity	0.052ng/mL
Entrez	110749
Expiration Date	Please enquire.